630-229

# ORTHOVISC® Sodium Hyaluronate FOR INTRA-ARTICULAR INJECTION

**ORTHOVISC®**  
**Sodium Hyaluronate**  
**FOR INTRA-ARTICULAR INJECTION**

**Non-surgical use**

**Sterile Injection**

**DESCRIPTION:**

ORTHOVISC® is a sterile, non-pyrogenic solution of sodium hyaluronate. ORTHOVISC® contains 15 mg/mL of sodium hyaluronate (NaHA) dissolved in physiological saline. The kinematic viscosity of the solution is adjusted to 20,000-70,000 centistokes, and its osmolality is approximately 340 milliosmoles.

**CHARACTERISTICS:**

Sodium Hyaluronate is a high molecular weight polysaccharide composed of sodium glucuronate and N-acetylglucosamine. Sodium hyaluronate is ubiquitously distributed throughout tissues of the body and is present in high concentrations in such tissues as vitreous humor, synovial fluid, umbilical cord and dermis.<sup>1</sup> Sodium hyaluronate functions as a tissue lubricant and is thought to play an important role in modulating the interactions between adjacent tissues.<sup>2,3</sup> It can also act as a viscoelastic support maintaining a separation between tissues. Different sodium hyaluronate preparations may have different molecular weights, but are thought to have the same chemical structure. The sodium hyaluronate in ORTHOVISC® has an average molecular weight greater than one million. ORTHOVISC® is non-inflammatory, is non-pyrogenic, and is well tolerated.<sup>4</sup> Sodium hyaluronate has been shown to be non-antigenic.<sup>4,5</sup> ORTHOVISC® does not interfere with normal wound healing processes.

**INDICATION:**

ORTHOVISC® is indicated for the symptomatic treatment of osteoarthritis of the knee.<sup>6</sup>

**DIRECTIONS FOR USE:**

ORTHOVISC® is administered by intra-articular injection one week apart for a total of three injections. Additional series of injections may be administered when clinically indicated. The required amount of ORTHOVISC® is slowly infused through a sterile hypodermic needle of suitable gauge into the selected joint space. A plastic-hubbed needle is recommended. The volume will vary depending upon the size of the joint space, not to exceed 2 mL for the knee.

**DO NOT OVERFILL JOINT SPACE.**

**CONTRAINDICATIONS:**

At the present time there are no known contraindications to the use of ORTHOVISC® when used in the symptomatic treatment of joint disease.

ORTHOVISC® may contain trace quantities of gram positive bacterial proteins and is contraindicated for patients with a history of such allergies. There are minimal risks associated with the procedure of injecting substances into joints in general; primarily infection and bleeding. Pre-existing infections of skin in the region of the intended injection or known systemic bleeding disorders may constitute relative or absolute contraindications.

**PRECAUTIONS:**

Precautions which are normally considered during injection of substances into joints are recommended. Only individuals familiar with accepted injection techniques for delivering agents to joint spaces should inject sodium hyaluronate for this application. An excess quantity of sodium hyaluronate is not to be used and the patient should be monitored closely. If pain increases during the injection procedure, the injection should be stopped and the needle withdrawn. The space should not be overfilled.

**ADVERSE REACTIONS:**

Sodium Hyaluronate is a natural component of the tissues of the body. ORTHOVISC® is thoroughly tested to determine that each batch is non-inflammatory. Adverse events reported in clinical trials have included mild or moderate reactions at the site of the injection. Based upon clinical studies, these reactions have occurred with a similar incidence in patients treated with saline injections. There is no evidence which attributes reactions to be specifically related to the product.

**HOW SUPPLIED:**

ORTHOVISC® is a sterile viscoelastic preparation supplied in a disposable glass syringe delivering 2 mL, appropriate for the knee, of sodium hyaluronate dissolved in physiological saline. Each mL of ORTHOVISC® contains 15 mg of sodium hyaluronate adjusted to 20,000-70,000 centistokes, 9 mg of sodium chloride and q.s. Sterile Water for Injection USP. ORTHOVISC® exhibits an osmolality of approximately 340 milliosmoles. Contents of unopened and undamaged syringe are sterile. Refrigerated ORTHOVISC® should be allowed to reach room temperature, approximately 20 to 45 minutes, prior to use.

**FOR INTRA-ARTICULAR USE. STORE BETWEEN 2°C AND 25°C. PROTECT FROM FREEZING.**

**CAUTION:** This product is available under medical recommendation and the administration must be performed by a physician.

**REFERENCES:**

1. Pruett, R.C., Schepens, C.L. and Swann, D.A. (1979) Hyaluronic Acid Vitreous Substitute. A Six Year Clinical Evaluation. Arch. Ophthalmol. 97, 2325.
2. Swan, D.A., Radin, E.L., Nazimiec, M., Weissner, P.A., Curran, N. and Lewinnek, G. (1974) Role of Hyaluronic Acid in Joint Lubrication. Ann. Rheum. Dis. 33, 318.
3. Radin, E.L., Paul, I.L., Swann, D.A. and Schottstadt, E.S. (1971) Lubrication of Synovial Membrane. Ann. Rheum. Dis. 30, 322.
4. Richter, W. (1974) Non-immunogenicity of Purified Hyaluronic Acid Preparations Tested by Passive Cutaneous Anaphylaxis. Int. Arc. Allergy 47, 211.
5. Richter, W., Ryde, E. M. and Zetterstrom, E.O. (1979) Non-immunogenicity of a Purified Sodium Hyaluronate Preparation in Man. Int. Arch. Appl. Immunol. 59, 45.
6. Peyron, J. (1993) Intra-articular Hyaluronan Injections in the Treatment of Osteoarthritis: State-of-the-Art Review. Journal of Rheumatology, Vol 20, Sup 39, 10-15.

ORTHOVISC® is a registered trademark of Anika Therapeutics, Inc. used under license by Pharmascience, Inc.

Manufactured by: Size  
Anika Therapeutics, Inc. 2 mL  
32 Wiggins Ave.  
Bedford, MA 01730 U.S.A.

Copyright© 2014 AML 500-269/B  
Anika Therapeutics, Inc. 03/2013

Imported by:  
PENDOPHARM, Division of Pharmascience Inc.  
Montreal (Quebec)  
Canada, H4P 2T4  
Medical Information: 1-888-550-6060, medinfo@pendopharm.com

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