REF M12D104T

**SportVis™** **EN**

#### INTENDED USES:

SportVis™ is intended to relieve pain and restore function in acute first- and second-degree ankle sprain and chronic (>3 months) lateral epicondylalgia (tennis elbow).

SportVis™ sodium hyaluronate (STABHA™: Soft Tissue Adapted Biocompatible Hyaluronic Acid) augments the sodium hyaluronate naturally present in the soft tissue surrounding damaged tendons and ligaments and provides support, lubrication, and hydration to the affected site. SportVis™ has been

clinically demonstrated to relieve pain and restore function in sports-related acute first- or second-degree ankle sprains and chronic lateral epicondylalgia (tennis elbow).

#### CONTRAINDICATION:

Patients with known sensitivity to sodium hyaluronate.

#### PRESENTATION:

SportVis™ is a clear solution of sterile 1% sodium hyaluronate in phosphate-buffered saline contained in a pre-filled syringe for peri-articular injection into the soft tissue surrounding tendons and ligaments.

Sodium hyaluronate is a

long-chain polysaccharide made up of repeating disaccharide units, which occurs naturally in the body.

SportVis™ has a pH and osmolality biocompatible with the soft tissue. 1.2 mL of SportVis™, sterilized by filtration, is enclosed within a disposable ready-to-use glass syringe, to be used with the appropriate gauge of needle. Needles are not included with SportVis™. Each syringe is contained within a blister pack. There are two blister packs in each box of SportVis™.

#### COMPOSITION:

SportVis™ consists of the following components: Water for injections, sodium hyaluronate (10 mg/mL),

sodium chloride, disodium hydrogen orthophosphate dihydrate, and sodium dihydrogen phosphate dihydrate.

#### DOSAGE AND ADMINISTRATION:

Peri-articular injection of SportVis™ should only be made by a healthcare professional trained in the specific peri-articular injection technique.

The contents of the syringe are sterile and should be injected using a sterile needle of an appropriate size for the injection site.

The area to be treated should be disinfected before injection.

Discard the syringe and needle after single use.

#### Please note the following injury-specific administration instructions:

##### Ankle sprains

Two injections, as follows:

- One 1.2 mL peri-articular injection of SportVis™, preferably within 48 hours of the first- or second-degree ankle sprain, **followed by a second injection** two to three days after the first injection.

- A **27-gauge** needle is recommended.

Peri-articular injections should be delivered during a single penetration along the anterior talofibular ligament using clinical landmarks. The injection is delivered in a fanning technique

along three planes from anteroposterior, medial, and lateral to the proximal ligamentous landmark. A third of the contents of each syringe should be released into the injury site during withdrawal of the syringe during each of the 3 phases.

##### Tennis elbow (lateral epicondylalgia)

Two injections, as follows:

- One 1.2 mL peri-articular injection of SportVis™ at the lateral elbow epicondyle site, **followed by a second injection** at the same site one week after the first injection.

- A **27-gauge** needle is recommended.

Identify the tenderest point of the epicondyle by gentle palpation. Position the

needle at 45 degrees to the point of maximal pain of the lateral epicondyle. After puncture of the skin, angle the needle parallel to the skin and insert it towards the point of maximal pain on the lateral epicondyle. The injection is to be administered in two planes using a fanning technique. Inject half the contents as the needle is withdrawn towards the skin **without exiting the skin.**

Rotate the needle 180 degrees (opposite direction) and insert the needle parallel to the skin towards the point of maximal pain on the lateral epicondyle. Inject the remaining contents as the needle is withdrawn towards the skin. Remove

the needle from the skin. Flex and extend the elbow five times and then internally and externally rotate five times.

#### WARNINGS AND PRECAUTIONS:

SportVis™ should only be injected by a healthcare professional trained in peri-articular injection techniques.

SportVis™ pre-filled syringes are single use only. The contents of the syringe should be used for one injection only. Any remaining sodium hyaluronate should be discarded. If a syringe is retained for a subsequent injection, there is a risk of contamination resulting in the possible infection of

the patient, and/or foreign body reaction.

SportVis™ should not be re-sterilized as the device performance may be compromised, which could cause serious harm to the patient's health and safety. SportVis™ must not be injected into blood vessels. Sodium hyaluronate has the potential to occlude the vessels, which could result in embolism or infarction. SportVis™ should not be injected into a hematoma.

Direct injection into tendons should be avoided, as this can lead to rupture. Use of an ultrasound probe to guide the injection will minimize this risk.

Do not inject into the soft tissue of patients if the

area of the injection site is infected or if there is evidence of acute or chronic skin disease in the area.

Sodium hyaluronate is manufactured by fermentation of *Streptococcus equi* and rigorously purified. However, the healthcare professional should consider the immunological and other potential risks that can be associated with the injection of any biological material. There is a risk of infection at the injection site, as with any peri-articular procedure.

There is no evidence of the safety of SportVis™ in patients who are pregnant or nursing.

The safety and effectiveness of SportVis™ have not been tested in children under 18 years of age.

Follow national or local guidelines for the safe use and disposal of needles. Obtain prompt medical attention if injury occurs. Do not use if the blister pack containing the syringe has been damaged. Do not use after the expiry date printed on the package.

#### ADVERSE REACTIONS:

Mild erythema that should resolve with time. General side effects that may occur during any periarticular injection include but are not limited to: blood

vessel occlusion, bleeding, infection, inflammation, nerve damage, pain, or device dislocation. Read the Warnings and Precautions and Dosage and Administration instructions before use.

#### INCOMPATIBILITIES:

SportVis™ has not been tested for compatibility with other substances for peri-articular injection. Therefore, the mixing or simultaneous administration with other peri-articular injectables is not recommended.

#### STORAGE:

Store between 5 °C and 25 °C. Do not freeze. Protect from light. Do not use if sterile packaging has been

damaged. Do not use after expiry date.

#### SYMBOLS USED ON THE PACKAGING:

	Lot number
	Expiry date
	For single use only
	Store between 5 °C and 25 °C
	Protect from light
	Do not freeze
	Do not use if sterile packaging has been damaged

 SportVis™ pre-filled syringes are aseptically filled with sterile filtered sodium hyaluronate

 Manufacturer
   
 Date of manufacture
   
 Full instructions, warnings, and precautions are included in this instruction for use leaflet.

 Catalogue number

**Pendopharm Medical Information:**
  
 1-888-550-6060
   
 medinfo@pendopharm.com
   
 Issue date: 11/2020

SportVis™ is a trademark of MDT Int'l. SA, used under licence.

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