INSTRUCTIONS FOR USE
MONOVISC®, Lightly Cross-Linked High Molecular Weight
Sodium Hyaluronate, 4 mL

DESCRIPTION:

MONOVISC® is a sterile, non-pyrogenic, translucent sodium hyaluronate solution, lightly cross-linked with a proprietary chemical cross-linker. MONOVISC® contains nominal 20 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline at physiological pH and osmolality. MONOVISC® is manufactured from ultra pure, high molecular weight sodium hyaluronate produced by bacterial fermentation. Hyaluronic acid is a natural complex polysaccharide of the glycosaminoglycan family.

CHARACTERISTICS:

Sodium Hyaluronate is a high molecular weight polysaccharide composed of sodium glucuronate and N-acetylglucosamine. The sodium hyaluronate in MONOVISC® is derived from bacterial fermentation. Hyaluronic acid is ubiquitously distributed throughout tissues of the body and is present in high concentrations in such tissues as vitreous humor, synovial fluid, umbilical cord and dermis. Sodium hyaluronate functions as a tissue lubricant and is thought to play an important role in modulating the interactions between adjacent tissues. It can also act as a viscoelastic support maintaining a separation between tissues. Different sodium hyaluronate preparations may have different molecular weights, but have the same chemical structure. The MONOVISC® cross-linked sodium hyaluronate injection is biocompatible, non-inflammatory and non-pyrogenic. Sodium hyaluronate preparations have been shown to be well tolerated in osteoarthritic synovial joints such as knee, hip, shoulder and ankle. The actions of MONOVISC are lubrication and mechanical support.

APPLICATION:

MONOVISC® Cross-Linked Sodium Hyaluronate Injection is a single, intra- articular injection of cross-linked sodium hyaluronate designed to treat the symptoms of osteoarthritis.

INDICATIONS:

MONOVISC® is indicated in the treatment of pain in osteoarthritis (OA) of the knee, hip, shoulder and ankle for patients who have failed to respond adequately to conservative non-pharmacologic therapy and to simple analgesics, e.g., acetaminophen.

DIRECTIONS FOR USE:

The required amount of MONOVISC® is injected through a sterile, disposable, hypodermic needle of suitable gauge into the selected joint space. The sterile needle should be attached to the MONOVISC® syringe by a health care professional using a health care facility-approved aseptic technique. Common needle gauges for injections into the knee are 18-21 gauge. The final needle selection for any procedure is determined by the physician. The health care provider should ensure proper penetration into joint synovial space prior to injecting MONOVISC®. The syringe is intended for single injection only; any unused product should be discarded per health care facility approved procedures.

DOSAGE:

MONOVISC® is a single dose preparation. The recommended dose is 4 mL (one syringe) per joint.

CONTRAINDICATIONS:

MONOVISC® is composed of cross-linked sodium hyaluronate and may contain trace amounts of gram positive bacterial proteins. The following pre-existing conditions may constitute relative or absolute contraindications to the use of MONOVISC®:

- Known sensitivity to hyaluronic acid based products
- Pre-existing infections of the skin region of the intended injection site
- Known infection of the index joint
- Known systemic bleeding disorders
- Known sensitivity to any of the materials contained in MONOVISC®

PRECAUTIONS:

- Precautions which are normally considered during injection of substances into joints are recommended.
- Only medical professionals trained in accepted injection techniques for delivering agents to joint spaces should inject sodium hyaluronate for this application.
- The amount of MONOVISC® necessary to be injected depends on specific site and patient anatomy and needs to be defined by the medical professional performing the procedure. An excess quantity of sodium hyaluronate should not to be used and the patient should be monitored closely.
- The synovial space should not be overfilled.
- If pain increased during the injection procedure, the injection should be stopped and the needle withdrawn.
- The safety and effectiveness of MONOVISC® concomitantly with other intra-articular injectables have not been tested.
- MONOVISC® has not been tested in pregnant or lactating women or in children.
- The effect of repeated injections of MONOVISC® has not been tested.

ADVERSE REACTIONS:

Hyaluronic acid is a natural component of the tissues of the body. MONOVISC® is thoroughly tested to determine that each batch conforms to the product quality attributes. MONOVISC® is proven to be non-inflammatory by biocompatibility testing.

Mild to moderate episodes of transient swelling and discomfort have occasionally been observed following intra-articular injection of sodium hyaluronate preparations. General risks associated with the procedure of injecting substances into joints may include infections and bleeding.


HOW SUPPLIED:

MONOVISC® is a sterile viscoelastic preparation supplied in a disposable glass syringe delivering 4 mL. Each mL of MONOVISC® contains 20 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline.

FOR INTRA-ARTICULAR USE. STORE BETWEEN 2°C AND 25°C.
PROTECT FROM FREEZING.

MONOVISC®, Lightly Cross-Linked High Molecular Weight
Sodium Hyaluronate, 4 mL

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AML 500-260 REV-C
2025-05

CAUTION: This product is available under medical recommendation and the administration must be performed by a physician.

MONOVISC® is a registered trademark of Anika Therapeutics, Inc. used under license by Pharmascience Inc.
Manufactured by:
Anika Therapeutics, Inc.
32 Wiggins Ave.
Bedford, MA 01730 U.S.A.

Imported by:
PENDOPHARM, Division of Pharmascience Inc.
Montreal (Quebec)
Canada, H4P 2T4
Medical Information: 1-888-550-6060, medinfo@pendopharm.com

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