**Eleviev™** Submucosal injectable composition 10 mL  
To facilitate endoscopic resection procedures  
10 mL per ampoule

#### INTENDED USES

Eleviev™ is a submucosal injectable composition intended for use in endoscopic resection procedures for submucosal lift of polyps, adenomas or other gastrointestinal mucosal lesions in the colon, prior to excision with a snare or endoscopic device.

Eleviev™ is injected into the submucosa, beneath the lesion to be excised. Eleviev™ is designed to elevate the mucosal layer and the tissue to be excised from the submucosal and the external muscular layer, providing a submucosal cushion of optimal height and duration, allowing the endoscopist an easier and safer resection procedure: endoscopic mucosal resection (EMR), hybrid EMR and endoscopic submucosal dissection (ESD).

Eleviev™ does not require any special apparatus or equipment, and it is designed to be used with the most common endoscopic resection devices.

The emulsion has physicochemical properties suitable for injection in the target site via a normal, commercially available endoscopic injection needle (e.g.: a 2.3 mm x 230 cm endoscopic injection needle) having a needle diameter of 23 gauge (23G) or less (not provided with Eleviev™).

The ampoules are supplied with a female Luer-lock closure to enable extraction of the emulsion with a common glass or plastic syringe with male Luer-lock connection fitting, to be subsequently connected to the endoscopic injection needle.

Following injection into the submucosa, Eleviev™ reconfigures to occupy the interstitial space and forms a colored submucosal cushion of optimal height that pushes the mucosa away from the submucosal layer, allowing for an easier resection procedure. Eleviev™, as a colored clear emulsion, helps to clarify the area where it is injected, to assist the endoscopist in visualizing the margins of the target lesion and performing the resection procedure, thereby decreasing the risk of damaging the external muscular layer, which could lead to perforation.

#### PRESENTATION

Eleviev™ is a sterile, colored clear emulsion.

Eleviev™ is sterilized by filtration and aseptically filled into a 10 mL polypropylene ampoule with a female Luer-lock closure that can easily be connected to a suitable disposable plastic or glass syringe with male Luer-lock connection fitting, to extract the emulsion and be injected through an endoscope via a normal, commercially available endoscopic injection needle (e.g.: a 2.3 mm x 230 cm endoscopic injection needle) having a needle diameter of 23 gauge (23G) or less (not provided with the device).

Eleviev™ is contained in 10 mL polypropylene ampoules pouched in an aluminum pouch under nitrogen, which is packed in a cardboard box.

#### COMPOSITION

The emulsion consists of the following components: Water for injection, Medium chain triglycerides, Poloxamer 188, Polyoxyl-15-Hydroxystearate, Sodium Chloride, Methylene blue.

#### DOSAGE AND ADMINISTRATION

Eleviev™ can be injected through an endoscope via a normal, commercially available endoscopic injection needle (e.g.: a 2.3 mm x 230 cm endoscopic injection needle) having a needle diameter of 23 gauge (23G) or less (not provided with Eleviev™).

The administered dose of Eleviev™ should be determined based on the dimensions of the lesion to be removed.

Inject into the submucosa the amount of Eleviev™ needed to form a submucosal cushion of optimal height and shape for the lesion to be removed.

During the procedure do not exceed a total dose of 50 mL per patient, either in single or multiple administrations.

#### PROCEDURE

- Open the aluminum pouch and detach one ampoule from the strip.
- Open the ampoule by twisting off the cap following the arrow direction.
- Screw a disposable, sterile plastic or glass syringe with male Luer-lock connection fitting (not provided with Eleviev™) directly to the ampoule.
- Turn the ampoule upside down, and extract the whole volume of the emulsion.
- Unscrew the disposable, sterile plastic or glass syringe from the ampoule, and expel the air from the syringe until some drops of the emulsion pour out of the tip.
- Securely screw the filled Luer-lock syringe to a suitable sterile endoscopic injection needle (not provided with Eleviev™).
- Fill the endoscopic injection needle with the emulsion and expel any air.
- Infuse a small quantity of the emulsion to determine the amount of force needed to inject the material.
- Retract the needle.
- Introduce the endoscopic injection needle through the working channel of the endoscope.
- Inject Eleviev™ immediately beneath the lesion to be excised in one or in a series of injections around the base of the lesion.
- Discard the empty ampoules. Store unused ampoules inside the aluminum pouch, inside the box, for up to 6 months from opening of the pouch.

#### CONTRAINDICATIONS

Patients with known sensitivity to any of the components contained in Eleviev™.

#### WARNINGS AND PRECAUTIONS

- The endoscopist injecting Eleviev™ must be experienced in the administration technique.
- The safety of Eleviev™ has not been established in pregnant or nursing women or in children under 18 years of age.
- Eleviev™ is provided in single-use ampoules. Eleviev™ should not be reused after first opening. Any emulsion not injected during the procedure should not be reused for another endoscopic procedure.
- Do not use if the primary packaging (ampoule) or secondary packaging (aluminum pouch) is damaged.
- Do not use if the twist-off cap is damaged.
- Do not use if the emulsion is not clear, shows any signs of opalescence or contains floating or precipitated visible particles.
- The product compatibility with other substances has not been tested.

#### ADVERSE REACTIONS

Rarely, local bleeding and/or inflammatory reaction could occur, which may or may not be associated with Eleviev™.

#### STORAGE

Store between 2°C and 25°C.  
Excursions permitted up to 40°C.  
Protect from light. Do not freeze.  
Do not use if packaging is damaged.  
Sterile product for single use only.

Do not use after expiry date.  
Revision n° 00, Date: April 2019

#### DESCRIPTION OF USED SYMBOLS

		
Manufacturer	Date of manufacture	Batch code
		
Catalogue number	Consult instructions for use	Do not reuse
		
Do not resterilize	Do not use if package is damaged	Sterilized using aseptic processing techniques

 Manufacturer:  
Cosmo Technologies Ltd.  
Riverside II, Sir John Rogerson's Quay  
Dublin 2, Ireland

Distributed by:  
Pendopharm, Division of Pharmascience Inc.  
Montreal, Canada, H4P 2T4

Medical Information: 1-888-550-6060  
medinfo@pendopharm.com  
www.pendopharm.com

™ Trademark of Cosmo Technologies Ltd.,  
used under license.

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Pendopharm, division de Pharmascience Inc.